



Valplast Care and Cleaning Instructions

To maintain maximum cleanliness and shine only **Val-Clean** solution should be used to soak and wash your Valplast in. Val-Clean can be purchased at Dream Dental or via ordervalplast.com.au

You will be provided with a container for storing and soaking your Valplast.

- Clean your Valplast appliance regularly. The Val-Clean Concentrated Denture Cleaner will safely, quickly, and effectively maintain cleanliness and shine
- Soak your Valplast appliance in Val-Clean for 10-15 minutes a day, or overnight at least three times a week
- To make up the Valplast solution you put $\frac{1}{4}$ of a tsp of solution to 150 ml of water
- Loose particles can be removed by placing the appliance under running water
- Keep your Valplast in water or in Val-clean whenever it is not being worn to keep the surface hydrated
- Remove your Valplast appliance prior to flossing and brushing your natural teeth

Please Do Not

- Do not use a toothbrush or toothpaste on a Valplast appliance, as this may remove the polish and roughen the surface over time
- Do not soak in anything other than Val-Clean
- Do not use other brands of denture cleaner as they are not effective and may contain ingredients that can damage the Valplast or its colour

Your Valplast must always be wet – If it's not in your mouth it must be kept in a container with Val-Clean or plain water

Should you experience any persistent irritation please contact **Dream Dental** on **3219 4020**